

2568 | PRT Scanner Module, 8-channel, 1000 Ω



产品概述: 2568 | PRT Scanner Module, 8-channel, 1000 Ω

For 1000-ohm PRTs, these modules provide all the same great features as the 2560 and 2562 Modules. The 2568 Module reads up to eight 1000-ohm PRTs and at 0°C is accurate to $\pm 0.01^\circ\text{C}$. Don't use an ohmmeter or multimeter to read your 1000-ohm PRTs when you can use a Black Stack loaded with convenient temperature functions.

型号



2568

PRT 扫描仪模块，8 通道，1000 Ω



Fluke. 让您的工作畅通无阻。

福禄克测试仪器（上海）有限公司 电话：400-810-3435 ©2025 福禄克公司
06/2025
北京福禄克世禄仪器维修和服务有限公司 电
话：400-615-1563 未经许可，本文档禁止修改
福禄克测试仪器（上海）有限公司上海维修中心 电
话：021-54402301, 021-54401908分机269
福禄克测试仪器（上海）有限公司深圳第一特约维修点
电话：0755-86337229